


# 물 품 규 격 서

설치장소	신공학관 3128호	수량	1
품명	고속분리 정제기		
	High Throughput Column Chromatography		
<p><b>세부내용 :</b></p> <p>1. 제조사 : Biotage</p> <p>2. 제품명 : ISO-1SV</p> <p>3. 주요사양 :</p> <ul style="list-style-type: none"> <li>* Specifications ;</li> <li>· Flow Rate : 1-200ml</li> <li>· Pressure limit : 145Psi (10bar)</li> <li>· Two Piston Gradient HPFC</li> <li>· Quatro Binary Gradient Pump</li> <li>· TLC-to-Step Gradient with Simulation</li> <li>· 3D Graphic for Complete Reveal</li> <li>· Co_Solvent Isocratic additive</li> <li>· Power : 100 - 240V, 50/60 Hz, 4.0 A</li> </ul> <p>Fraction Collector</p> <ul style="list-style-type: none"> <li>· Manual &amp; Automatic Moving Arm</li> <li>· Standard : Fraction Volume : Max 4.8 L</li> <li>· 16mm x 150mm (4ea) 분리형 독립 RACK</li> </ul> <p>Fraction Modes :</p> <ul style="list-style-type: none"> <li>· Volume, Threshold, Threshold with Volume</li> <li>· Low Slope, Medium Slope, Custom Slope</li> <li>· Real-Time PDA Scanning (200-400nm)</li> <li>· UV Wavelength (with BaseLine Correction)</li> <li>· 2D and 3D Spectral Analysis (λ- All Detection)</li> </ul> <p>4. 사진 :</p> <div style="text-align: center; margin: 20px 0;">  </div>			
<p><b>용도설명:</b> 자동 컬럼크로마토그래피 기기로 혼합물을 알맞는 전개용매 조건에서 자동으로 분리</p>			
<p><b>기타참고자료 :</b> 업체 카달로그</p>			